


Mirror Chart

MC-5/MC-5S



MC-5/MC-5S

MC-5/MC-5S supports you in the path to adapting quickly to new social distance recommendations whilst providing sophisticated visual acuity (VA) and space saving.



Installation distance: **40cm***1



Installation depth: **21cm**




MC-5S

NEW refraction system with "Circular polarization"


Circular polarization

The built-in circular polarization panel provides stability even if an individual moves their head. This is due to the stability provided by the binocular balance visual function test with a polarization chart.*2*3

3D contents can also be performed on the MC-5S.



OK!



OK!

MC-5

NEW LCD Remote controller !

LR-5

The new LCD remote control is more comfortable to hold due its new built-in LR-5 alignment emitter.



Alignment emitter



Points!

The infrared transmitter is now wider*4.

*1 Minimum installation distance 40cm with MC-5, 60cm with MC-5S.

*2 With wearing the circular polarization glasses while the test.

*3 The Circular polarization glasses are not a standard accessory while wearing circular polarization glasses.

*4 Compared with LR-4.


VA test meets social distance recommendations

Remote Tablet operation enables social distance.^{*5}

*5 in case MC-5 is connected with CV-5000PRO and the tablet PC is utilized as the operation controller.

Space saving VA test and flexible installation

- The VA test with a small installation distance.
- Flexible placement positioning in your practice or clinic.
- The multiple refraction lanes in a small area.



Optional accessory for MC-5



LCD remote controller
LR-5

Combination devices with MC-5S



Compu Vision
CV-5000PRO



1 dial controller
KB-50S

Specifications

Installation distance	MC-5 : 40cm, 65cm, 85cm, 100cm, 105cm, 155cm MC-5S : 60cm, 80cm, 95cm, 110cm, 150cm*	*Only for H1 type
VA value	MC-5 : 0.03 ~ 2.0** MC-5S : 0.03 ~ 2.0	**0.03 can be displayed when MC-5 is connected with CV-5000PRO.
Test chart	MC-5 : Landolt ring, numerals tumbling E, R&G, astigmatism, etc. MC-5S : Landolt ring, numerals, alphabets, tumbling E, R&G, astigmatism, etc.	
Dimensions	450mm (W), 207mm (D), 621mm (H)	
Weight	MC-5 : 13.6kg, MC-5S : 13.1kg	
Power supply	AC: 100V-240V, 50-60Hz	

Chart sample

(F type for MC-5)

												4 5 8 2 9 3	5 4 3 6 7 6 8 2 9 4 3 2 8 6 5	2 8 5 6 3 9 3 4 2 7 5 6 9 8 4	9 4 5 2 7 6 8 3 4 9 5 9 2 8 3							

IMPORTANT In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.
Not all products, services, or offers are available in all markets. Contact your local distributor for country-specific information and availability.

TOPCON INSTRUMENTS (MALAYSIA) SDN. BHD.
(Regional Office for Topcon Healthcare Southeast Asia)
Unit 2, 4, Jalan Pensyarah U1/28, Hicom-glenmarie Industrial Park,
40150 Shah Alam, Selangor, MALAYSIA
Phone: +603-766 16260 Fax: +603-766 16261
Email: mys_tim_marketing_sm@topcon.com
www.topconhealthcare.my

TOPCON SINGAPORE MEDICAL PTE. LTD.
100G Pasir Panjang Road, #02-18, Interlocal Centre,
SINGAPORE 118523
Phone: +65-68720606 Fax: +65-67736150
E-mail: med.sales.sg@topcon.com
www.topconhealthcare.sg

TOPCON INSTRUMENTS (THAILAND) CO., LTD.
77/162 Sinnsathorn Tower, 37th Floor, Krungthoburi Rd.,
Klongtomsai, Klongsarn, Bangkok 10600, THAILAND
Phone: +66-02-440-1152 Fax: +66-02-440-1158
Email: tha_medical@topcon.com
www.eyecare.topcon.co.th

MEHRA EYETECH PRIVATE LIMITED
801 B Wing, Lotus Corporate Park, Graham Firth Steel Compound
Goregaon (East) Mumbai 400063 Maharashtra, INDIA
Phone: +91-22-61285455
E-mail: sales@mehraeyetech.in
www.topconhealthcare.in

TOPCON (BEIJING) MEDICAL TECHNOLOGY CO., LTD.
Room 2808, Tower C, JinChangAn Building, No.82, Middle Section of
East 4th Ring Road, Chaoyang District, Beijing 100124, P.R. CHINA
Phone: +86-10-8794-5176
E-mail: cn_marketing@topcon.com
www.topcon-china.net/

TOPCON KOREA MEDICAL CO., LTD.
2F YK Building, 205, Dogok-ro, Gangnam-gu,
Seoul, Republic of Korea
Phone: +82-2-6959-7947
E-mail: tkm@topcon.com
www.topconhealthcare.kr/



TOPCON CORPORATION

75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, JAPAN.
Phone: +81-(0)3-3558-2522/2502
Fax: +81-(0)3-3965-6898
www.topconhealthcare.jp

